

Rat Probe 4.0 Data Acquisition Pseudo code

	<pre>//Setup 1) Include files 2) Configure Bluetooth 3) Configure ADCs 4) Configure other required peripherals 5) Global variable declarations 6) Function prototypes 7) Other required initialization 8) Initialize all LEDs = off and scope trigger = low 1) int PD1_Ch1_730nm // PD1 measured with 730nm LED on channel 1 illuminated 2) int PD2_Ch1_730nm // PD2 measured with 730nm LED on channel 1 illuminated 3) int PD1_Ch1_850nm // PD1 measured with 850nm LED on channel 1 illuminated 4) int PD2_Ch1_850nm // etc. 5) int PD2_Ch2_730nm // etc. 6) int PD1_Ch2_730nm // etc. 7) int PD2_Ch2_850nm // etc. 8) int PD1_Ch2_850nm // etc.</pre>
	<pre>// Make the Measurements 9) Loop forever 10) i = 0 11) PD1_Ch1_730nm = 0 //Initialize the photodiode measurement variables 12) PD2_Ch1_730nm = 0 13) PD1_Ch1_850nm = 0 14) PD2_Ch1_850nm = 0 15) PD2_Ch2_730nm = 0 16) PD1_Ch2_730nm = 0 17) PD2_Ch2_850nm = 0 18) PD1_Ch2_850nm = 0</pre>

19) Loop as many times as possible in < 1s

20) Scope trigger pin = high // Trigger scope

21) Scope trigger pin = low

22) Call Read_Photodiodes(Ch1_730nm_LED, PD1, PD2)

23) PD1_Ch1_730nm += PD1

24) PD2_Ch1_730nm += PD2

25) Call Read_Photodiodes(Ch1_850nm_LED, PD1, PD2)

26) PD1_Ch1_850nm += PD1

27) PD2_Ch1_850nm += PD2

28) Call Read_Photodiodes(Ch2_730nm_LED, PD1, PD2)

29) PD2_Ch2_730nm += PD2

30) PD1_Ch2_730nm += PD1

31) Call Read_Photodiodes(Ch2_850nm_LED, PD1, PD2)

32) PD2_Ch2_850nm += PD2

33) PD1_Ch2_850nm += PD1

34) i++

35) End Loop

36) PD1_Ch1_730nm = PD1_Ch1_730nm/i // Calculate averages

37) PD2_Ch1_730nm = PD2_Ch1_730nm/i

38) PD1_Ch1_850nm = PD1_Ch1_850nm/i

39) PD2_Ch1_850nm = PD2_Ch1_850nm/i

40) PD1_Ch2_730nm = PD1_Ch2_730nm/i

41) PD2_Ch2_730nm = PD2_Ch2_730nm/i

42) PD1_Ch2_850nm = PD1_Ch2_850nm/i

43) PD2_Ch1_850nm = PD2_Ch1_850nm/i

44) Transmit data to laptop via Bluetooth

45) Delay such that the outer loop executes in precisely 1s.

46) End Loop

	<pre> // // Read_Photodiodes Subroutine (24ms execution time) // Read_Photodiodes(Selected_LED <i>input parameter</i>, P_Diode_1 <i>output parameter</i>, P_Diode_2 <i>output parameter</i>) // Read the Dark Signals Before Making any Measurement 1) int D_Pdiode_1, D_Pdiode_2 // Dark PD measurement variables 2) int I_Pdiode_1, I_Pdiode_2 // Illuminated PD measurement variables 3) All LEDs = off 4) D_Pdiode_1 = 0, D_Pdiode_2 = 0 5) k = 0 6) Delay for 10ms //Ensure the preamps have settled to 0.25V 7) 8) Loop as many times as possible in < 2ms 9) D_Pdiode_1 += New ADC reading on channel 1 10) D_Pdiode_2 += New ADC reading on channel 2 11) k++ 12) End 13) D_Pdiode_1 = D_Pdiode_1/k // Average of Dark Photodiode 1 14) D_Pdiode_2 = D_Pdiode_2/k // Average of Dark Photodiode 2 // Read the illuminated Signals 15) Selected_LED = on 16) I_Pdiode_1 = 0, I_Pdiode_2 = 0, k = 0 17) Delay for 10ms // Allow preamps to settle 18) Loop as many times as possible in < 2ms 19) I_Pdiode_1 += New ADC reading on channel 1 20) I_Pdiode_2 += New ADC reading on channel 2 21) k = k + 1 22) End 23) Selected_LED = off 24) I_Pdiode_1 = I_Pdiode_1/k // Average of illuminated Photodiode 1 25) I_Pdiode_2 = I_Pdiode_2/k // Average of illuminated Photodiode 2 26) P_Diode_1 = I_Pdiode_1 - D_Pdiode_1 // Subtract dark offset 27) P_Diode_2 = I_Pdiode_2 - D_Pdiode_2 // Subtract dark offset 28) Return (P_Diode_1, P_Diode_2) </pre>
--	--